

DENTAL TRIBUNE

—The World's Dental Newspaper • Middle East & Africa Edition—

PRINTED IN DUBAI

www.dental-tribune.me

JANUARY-FEBRUARY 2014 | NO. 1, VOL. 3

EVENTS

"The Saudi Dental Society's main goal is ..."

>SDS IDC 2014 Insertion



INDUSTRY

AEEDC Dubai
04-06.02.2014

>Page 38



ACADEMIA TRIBUNE

College of Dentistry -
Sharjah University -
welcomes new dean

>Page 32



Ten years of Dental Tribune International

By Dental Tribune International

LEIPZIG, Germany: On Monday, 9 December, Dental Tribune International (DTI) celebrated another milestone in its history in the dental publishing business. The date marked the tenth anniversary of the company's foundation in 2003.

"What started about 20 years ago in Germany has developed into an international dental

and educational network over the past decade. Although the very first Dental Tribune edition was published in 1893, the real globalisation of the business ultimately started with the involvement of our publishing partners around the world. Today, we are able to produce unique local content with a truly global perspective," said DTI's CEO Torsten Oemus.

Currently, DTI offers more than 130 print publications and multiple websites that reach over



Torsten Oemus Publishers' Meeting at IDS Cologne 2013 (Credit: Daniel Zimmermann, DTI)

650,000 dentists in more than 90 countries in 27 languages. Since its foundation in 2003, the publisher has become the official media partner of a number of major events dedicated to the dental industry, such as the International Dental Show, the

Greater New York Dental Meeting, the annual congress of the FDI World Dental Federation and IDEM Singapore.

In addition to its print and online publications, the group can look back on its projects

in dental education and its flagship e-learning platform, Dental Tribune Study Club, in particular. Since its foundation in 2009, the sophisticated Web-

> Page 2

Meeting review: 89th Greater New York Dental Meeting



GNYDM in New York (Photo Courtesy of DT America)

By Dental Tribune America

NEW YORK, USA: The 89th Greater New York Dental Meeting was held Nov. 29 to Dec. 4 at the Jacob K. Javits Center in Manhattan, offering meeting attendees the opportu-

nity to expand their professional knowledge and expertise, and to visit with hundreds of exhibiting companies to learn about new products and services. Educational sessions covered all the bases — and all the body parts.

All-day live sessions at the Live Dentistry Arena included Dr. Jack Griffin Jr. demonstrating dependable, efficient preparations for monolithic lithium disilicate or zirconia crowns, digital impressions, cementation and finishing, featuring some of the most dependable materials available today. Dr. Aeklavya Panjali also spoke in a Live Dentistry session, as he surgically placed and restored a complex immediate placement implant case.

Orasoptic's Vanessa Velasco and Tom Lindsey brought the best in human engineering to the Dental Tribune Media Lounge with a presentation on two of the company's most recent additions: the groundbreaking XVI all-in-one dental loupe and headlight and the Body Guard PRO saddle chair. Velasco and Lindsey also talked about the Body Guard Pro saddle seat, which they described as stunning and comfortable.

At the Laser Pavilion Lecture Series, Dr. William R. Gianni of

> Page 40

UPTO DATE

PREVENTION KNOWLEDGE CONTINUED EDUCATION

Scientific Exchange Seminar for the whole dental team

The Role of the Dental Hygienist in Comprehensive Non-Surgical Periodontics

Prof. Crawford Bain

19:00 - 22:30
Friday 31st January 2014
Gloria Hotel, Dubai, UAE

Free Registration at
dtmea@dental-tribune.me

2 CME HOURS

ADA CERP
CAPP designates this activity for 2 continuing education credits.

New hyaluronic matrix accelerates soft-tissue healing



The new Hyaloss matrix. (Photo courtesy of imperiOs)

By Dental Tribune International

FRANKFURT, Germany: German specialist in the field of bone augmentation imperiOs has presented a

new product that is based on hyaluronic acid, a substance primarily used in plastic surgery but with potential applications in dentistry. According to the company, Hyaloss matrix, a hyaluronic matrix, promotes and accelerates the healing process to a significant degree.

According to imperiOs, Hyaloss matrix is a bioactive and resorbable matrix composed of hyaluronic acid fibres and is produced through esterification of the hyaluronic acid molecule with benzyl alcohol. Once the fibres of the matrix come into contact with liquid, the matrix gelatinises and can easily be inserted into

the respective bone defects, where it releases hyaluronic acid gradually.

Through activation of angiogenesis and mesenchymal stem cells, the matrix promotes regeneration processes during the first ten days after surgery in particular and thus contributes to faster healing.

According to the company, Hyaloss matrix is recommended for use in intraosseous and periodontal defects. The best results can be achieved when the matrix is mixed with autologous bone grafts, imperiOs stated.

An advantage of the matrix is that it can be stored at room temperature and can thus be used immediately to fill defects. In smaller periodontal defects, it can even be used unmixed. [DTI](#)

> Page 1

based continuing-education portal has become an essential part of DTI's product portfolio and serves more than 180,000 members worldwide. Today, DTI also provides continuing medical education through its clinical master's programme, Tribune CME, which offers comprehensive training in aesthetics, orthodontics and implantology, among other fields of dentistry.

With its headquarters based in Leipzig, Germany, the publisher today has a total of 34 licence

partners covering over 90 dental markets, including China, Brazil, South Africa, the US and the Middle East. Only recently, a new partner in Israel joined the publishing group. Talks for future projects in Sweden and Ukraine are currently underway.

The group is looking forward to celebrating the anniversary with its partners at its Annual Publishers' Meeting, which will be held in Turin in July of next year. [DTI](#)



Publishers' Meeting at IDS Cologne 2013
(Credit: Daniel Zimmermann, DTI)

REDUCES AND HELPS PREVENT GUM PROBLEMS IN 4 WEEKS TO ESCAPE THE CYCLE OF GINGIVITIS



Crest

Oral-B

CLINIC LINE
GUM PROTECTION SYSTEM

CLINICALLY PROVEN TO REDUCE GUM PROBLEMS

continuing the care that starts in your chair



International Imprint

Licensing by Dental Tribune International Group

Editor	Daniel Zimmermann newsroom@dental-tribune.com +49 341 48 474 107
Clinical Editor	Magda Wojtkiewicz
Online Editors	Yvonne Bachmann Claudia Duschek
Copy Editors	Sabrina Raaff Hans Motschmann
Publisher/President/CEO	Torsten Oemus
Director of Finance & Controlling	Dan Wunderlich
Business Development Manager	Claudia Salwiczek
Media Sales Managers	Matthias Diessner (Key Accounts) Melissa Brown (International) Peter Witteczek (Asia Pacific) Maria Kaiser (USA) Weridiana Mageswki (Latin America) Hélène Carpentier (Europe)
Marketing & Sales Services	Esther Wodarski
Accounting	Karen Hamatschek
Executive Producer	Gernot Meyer

Dental Tribune International
Holbeinstr. 29, 04229 Leipzig, Germany
Tel: +49 341 4 84 74 302 |
Fax: +49 341 4 84 74 173
www.dental-tribune.com
info@dental-tribune.com

Regional Offices

Asia Pacific
Dental Tribune Asia Pacific Limited
Room A, 20/F, Harvard Commercial Building,
105-111 Thomson Road, Wanchai, Hong Kong
Tel: +852 3113 6177 | Fax: +852 3113 6199

The Americas

Tribune America, LLC
116 West 23rd Street, Ste. 500, New York,
N.Y. 10011, USA
Tel: +1 212 244 7181 | Fax: +1 212 224 7185

Dental Tribune Middle East & Africa

Editorial Board

Dr. Aisha Sultan Alsuwaidi, President Emirates Dental Society, UAE
Dr. Nabeel Humood Alsabeeha, Consultant Prosthodontist,
Ministry of Health, UAE
Dr. Dobrina Mollova, MDS, UAE
Dr. Munir Silwadi, BDS, MRCDSO, DUSS, FADI, FICD, UAE
Dr. Khaled Abouseada, BDS, MS, orthod. cert.
Consultant Orthodontist, KSA
Dr. Rabih Abi Nader, B.D.S, D.E.S, D.E.S.S, Oral Surgeon and Implantologist, UAE
Rodny Abdallah, Dental Technician, CDT, Lebanon

Partners

Emirates Dental Society
Saudi Dental Society
Lebanese Dental Society

Director of mCME

Dr. Dobrina Mollova
info@cappmea.com
+971 50 42 43072

Business Development Manager

Tzvetan Deyanov
deyanov@dental-tribune.com
+971 55 11 28 581

Elegant design and easy handling are a winning combination.



CEREC OMNICAM

THE EVOLUTION OF SIMPLICITY



The new CEREC Omnicam combines powder-free ease of handling and natural color reproduction to provide an inspiring treatment experience. Discover the new simplicity of digital dentistry – exemplified by Sirona's premium camera portfolio: CEREC Omnicam and CEREC Bluecam. **Enjoy every day. With Sirona.**

UNRIVALLED HANDLING ■ POWDER-FREE ■ SCANNING IN NATURAL COLOR

Dubai Contact: SIRONA Dental Systems Ltd · Dr Amro Adel · Building 49, Suite 304
Dubai Healthcare City · Telephone: + 971 4 375 2355 · E-Mail: amro.adel@sirona.com

Kuwait Contact: SIRONA Dental Systems GmbH · Dr. Mostafa Al.Khouly
Office: + 965 2 224 6063 · Mobile: + 965 9 800 2225

sirona.com



The Dental Company

sirona.

“Dubai School of Dental Medicine has been highly successful”



By Dental Tribune Middle East

DUBAI, UAE: A year after the successful opening of the Dubai School of Dental Medicine we managed to catch up with Professor David Wray, Dean of DSDM to learn more about the developments since the opening.

DTMEA: Professor Wray, almost a year ago you opened the doors of the Dubai School of Dental Medicine with the aim to attract dental students and provide them with the oppor-

tunity to “Embark on a three year intensive clinical training program with a clear didactic component and a research dissertation”. How has the school developed since the opening?

Professor Wray: Since opening to students in January of this year DSDM has been highly successful. We accepted residents into our Pediatric Dentistry and Orthodontic programs in January and after completing our initial accreditation process with the CAA we have now accepted residents into all the specialty programs which we offer including Oral Surgery, Periodontology, Prosthodontics and Endodontics.

Seeing the growth of the DSDM since its launch, what sets you apart from the other Dental Universities in the UAE and the region?

The big bonus about DSDM is that graduates, as well as receiving their Master's Degree at the end of their course, also receive a conjoint Specialty

Membership Diploma from the Royal College of Surgeons of Edinburgh. This Specialty Membership Diploma gives

“DSDM is part of Dubai Healthcare City”

DSDM graduates a clinical qualification to complement their academic degree and a membership to the Royal College which is globally recognized, and DSDM provides a stimulating environment for graduates to carry out their postgraduate degree.

What are the major reasons for the regional dental students to choose DSDM?

DSDM provides both academic and clinical training within Dubai without having to leave the Gulf region which is a huge benefit especially to female residents with domestic commitments.

This can be achieved without

compromise since we have a world-class international faculty of teachers. DSDM provides education and training of the highest quality in a range of specializations to its students, and has also received initial academic accreditation from the Commission for Academic Accreditation of the UAE Ministry of Higher Education and Scientific Research. The school is led by world-renowned specialists and is home to some of the world's top specialists in the field. In addition, DSDM's Partnership with The Royal College of Surgeons of Edinburgh (UK) provides postgraduate students, through a conjoint exam, with a membership diploma which guarantees recognition internationally at specialty level in all areas. DSDM is part of Dubai Healthcare City, which is the world's first healthcare free-zone, and DSDM students are able to make the most of the world-class education facilities offered by the medical hub.

What are the further plans of DSDM in the coming 2 to 3 years?

DSDM has plans to expand its clinical facilities to cope with the increased numbers of residents and patients expected next year. Our partnership with the Royal College of Surgeons of Edinburgh is also progressing and we have already run four college examinations in DHCC this year. We anticipate that we will become one of the largest global hubs for dental postgraduate examinations in the near future.

Could you please share your thoughts on the level of Dentistry in the UAE and the region?

The quality of dental education and clinical practice in the UAE and the region is already very high but we hope to drive standards even higher with our postgraduate programs and of course, our graduates will be benchmarked clinically with the global standards set by the Edinburgh Royal College.

“We have now accepted residents into all the specialty programs”

Could you emphasize on the dental industry developments in the region?

We have seen dental companies expanding within the region and new ones coming. This reflects the interest in serving the dental community and the needs of the patients. In addition, the many conferences which Dubai hosts act as a platform to showcase the latest technology and how it can benefit the patients.

Through your experience as an educator, what should young students target when becoming

a dental professional?

Dentistry is largely a community-based profession and we need both general dentists and specialists to provide a full range of healthcare support for our patients so there is not just one path for a young dentist to follow. DSDM is dedicated to training specialists but there is a hugely important role for the generalist. Regardless of the path a young dentist chooses, it is important not to lose sight of the need for lifelong learning through continuing professional development be it in formal courses or through attendance at courses and self-directed learning. Dentists have a pro-

“DSDM provides both academic and clinical training within Dubai”

fessional responsibility to keep up-to-date and DSDM is here to support them in that endeavor.

How can you describe the dental students in Dubai?

The students in Dubai today are wonderful, mature, professional young dentists who are a family working towards a common goal. Our residents are a joy to teach and a privilege to be colleagues with.

Do you have any tips you would like to share with the young dental professionals?

Because of the explosion of new knowledge, young dentists must be constantly diligent and always treat the patients holistically. We don't just do fillings; we care for the total oral health needs and wellbeing of our patients. We should be proud of our profession and the care we can provide.

Do you have anything else to share with the dental audience?

Because DSDM is dedicated to specialization in dentistry we now represent a fantastic facility providing secondary oral healthcare to the community. Many dentists in UAE do not have the facilities or capacity to provide a fully comprehensive level of care and DSDM is here

to provide support and expertise to the dental community.

Whether the patients have complex restorative problems or serious mucosal disease, DSDM is here to help and welcomes referrals from all branches of the profession. **DT**

KaVo DIAGNOcam

Images that change your world.

KaVo DIAGNOcam – a whole new perspective on caries

- Significantly improved diagnosis quality – in unsurpassed image quality
- Ideal for patient information and outstanding monitoring
- X-ray-free imaging method for caries identification

KaVo DIAGNOcam – simply illuminating

Find out more about KaVo DIAGNOcam: www.kavo.de/diagnocam

KaVo. Dental Excellence.

KaVo Dental GmbH · Arjaan Tower 9th Floor · Dubai Media City, UAE · PO Box 71569
Phone: +971 4 433 21 86 · Fax: +971 4 457 93 73 · Email: info.mea@kavo.com · www.kavo.com/MEA

Contact Information

Direct: +971 4 424 8705
Fax: +971 4 424 8686
Dubai Healthcare City,
Bldg 34
Dubai, UAE
www.dsdm.ae



دبي مدرسة
لطب الأسنان
DUBAI SCHOOL
OF DENTAL MEDICINE



تقدم عيادة دبي لطب الأسنان عروض متكاملة لمعالجة وتنظيف الأسنان، حيث سيحرص نخبة من أطبائنا المحترفين على إضافة التألق لابتسامتك.

اختصاصات عيادة دبي لطب الأسنان:

تبييض الأسنان | تقويم الأسنان | علاج أمراض اللثة | زراعة الأسنان | التعويضات السنية
حشوات السيراميك | طب أسنان الأطفال | علاج الجذور | جراحة الفم

Dubai Dental Clinic offers treatments across all aspects of regular dental and hygienist services.
Our team of experts bring their experience to create a sparkle in your smile.

Dubai Dental Clinic specializes in:

Teeth Whitening | Orthodontics | Periodontal Treatments | Veneers | Dental Implants
Crowns | Ceramic Fillings | Pediatric Dentistry | Root Canals | Oral Surgery

لمزيد من المعلومات يرجى الاتصال على ٨٠٠ - DENTAL أو زيارة موقعنا الإلكتروني
Call us on 800 - DENTAL or visit our website for more information

الهاتف ٨٧٠٣ ٤٢٤ ٩٧١ | الفاكس ٨٦٨٦ ٤٢٤ ٩٧١ | البريد الإلكتروني clinic@dsdm.ac.ae
مدينة دبي الطبية، مبنى ٣٤، الطابق الأرضي
تفتح العيادة من السبت إلى الخميس من الساعة ٩ صباحاً وحتى الساعة ٦ مساءً

Phone +971 4 424 8703 | Fax +971 4 424 8686 | E-mail: clinic@dsdm.ac.ae
Dubai Healthcare City, Building 34, Ground Floor
Clinic hours Saturday to Wednesday from 9:00 a.m. to 6:00 p.m.
www.dsdm.ac.ae



كلية دبي لطب الأسنان
DUBAI SCHOOL
OF DENTAL MEDICINE

Mineral trioxide aggregate revisited: a cement for all seasons



mCME articles in Dental Tribune have been approved by HAAD as having educational content for CME credit hours. This article has been approved for 2 CME credit hours.

ADA C-E-R-P® Continuing Education Recognition Program

Centre for Advanced Professional Practices (CAPP) is an ADA CERP Recognized Provider. ADA CERP is a service of the American Dental Association to assist dental professionals in identifying quality providers of continuing dental education. ADA CERP does not approve or endorse individual courses or instructors, nor does it imply acceptance of credit hours by boards of dentistry.

CAPP designates this activity for 2 CE credits.

By Gary Glassman, DDS, FRCD

Pulpal and periradicular pathology develop when the dental pulp and periradicular tissues become exposed to microorganisms. In experimental, germ-free conditions, pulpal and periradicular tissues fail to show the development of pathosis and associated lesions when exposed to bacteria.^{1,2} The conclusion: Microorganisms are the main irritants of the dental pulp and periodontium, and sealing the pathways of communication between the root canal system and the periradicular tissues is imperative if bacterial leakage is to be prevented.

An ideal orthograde or retrograde filling material that seals the pathways of communication between the root canal system and its surrounding tissues should be non-toxic, non-carcinogenic, biocompatible, insoluble in tissue fluids and dimensionally stable.^{3,4} Furthermore, the presence of moisture should not affect its sealing ability; it should be easy to use and be radiopaque for recognition on radiographs.⁴

Because existing restorative materials used in endodontics did not possess these “ideal” characteristics, 4 mineral trioxide aggregate (MTA) was developed and recommended initially as a root-end filling material and subsequently has been used for pulp capping, pulpotomy, apexogenesis, apical barrier formation in teeth with open apices, repair of root perforations and, most recently, in revascularization cases. MTA has been recognized as a bioactive material.^{5,6}

MTA has been shown to seal off the pathways of communication between the root canal system and surrounding tissues, significantly reducing bacterial migration.⁷ It is made up of fine hydrophilic particles that set in the presence of water, and it is composed of tricalcium silicate, dicalcium silicate, tricalcium aluminate, tetracalcium aluminoferrite, calcium sulfate dihydrate (gypsum) and bismuth oxide, which provides it with radiopacity.⁸

Portland cement is the most common type of cement in general use around the world, used as a basic ingredient of concrete, mortar, stucco and most nonspecialty grout. It usually originates from limestone.

MTA is available as gray MTA and white MTA. The crystalline structure and chemical composition of gray and white MTA are similar, except for the presence of iron in gray MTA. Both contain bismuth oxide and calcium silicate oxide. Portland cement is composed mainly of calcium silicate oxide and does not contain bismuth oxide but does contain potassium. Calcium oxide is added in both Angelus white and gray MTA (Angelus, Londrina, Brazil) to reduce the setting time, which is too long in MTA cements of other brands (Fig. 1).

MTA has a similar mechanism of action to calcium hydroxide⁹ in that the main component of the material, calcium oxide, when in contact with a humid environment, is converted into calcium hydroxide.¹⁰ This results in a high pH of 12.5, making its surroundings inhospitable for bacterial growth and producing an antibacterial effect for a long period of time. But unlike calcium hydroxide products, such as Dycal® (DENTSPLY, York, Pa.) and MTA Angelus (Angelus, Londrina, Brazil), it has very low solubility, so it maintains a hard, excellent marginal seal.

Finally, unlike most dental materials, MTA actually needs moisture to set, so it thrives in a moist environment. Of the commercially available MTA products, MTA Angelus is well suited for most of the indicated endodontic procedures due to its setting time of 10 minutes, compared with the four-hour setting time of the other commercially available MTA. It is also packaged in air-tight bottles, allowing the practitioner to use only what is exactly needed, without introducing undue moisture into the remainder and without waste.¹¹

Endodontic revascularization

Treatment of the immature, non-vital tooth with apical pathology presents several challenges. The mechanical cleaning and shaping of such a tooth with a blunderbuss canal is difficult, if not impossible, to achieve predictably. The thin, fragile lateral dentinal walls can fracture during mechanical filing, and the large volume of necrotic debris contained in a wide root canal is difficult to completely disinfect.¹²

A new technique is presented to revascularize immature permanent teeth with apical

periodontitis. The canal is disinfected with copious irrigation and a combination of three antibiotics. After the disinfection protocol is complete, the apex is mechanically irritated to initiate bleeding into the canal to produce a blood clot to the level of the cemento-enamel junction.

A double seal of the coronal access is then made, first with MTA over the blood clot and then a bonded composite. The combination of a disinfected canal, a matrix into which new tissue could grow, and an effective coronal seal appears to have the ability to produce an environment necessary for successful revasculariza-

A case of mistaken identity

A 15-year-old girl of Asian descent was referred to the author's private endodontic clinic for evaluation on the lower left second premolar. The healthy young patient with an unremarkable medical history presented with a history of buccal swelling of the left mandibular area and discomfort to direct pressure on the tooth. On clinical examination, the patient was asymptomatic, and the tooth appeared intact, without caries. The presence of an enamel pearl on tooth #45 suggested that one may have been present on this tooth, which was frac-

testing, with mild sensitivity on percussion and palpation. Because of the presence of a wider than 4 mm open apex and thin dentinal walls prone to possible future fracture,¹⁴ it was felt that an attempt to achieve regeneration of the pulp should be made by a technique similar to that described by Rule and Winter¹⁵ and Iwaya et al.¹⁶

An access cavity was made, purulent hemorrhagic drainage obtained, and the necrotic nature of the pulp confirmed. The root canal was slowly flushed with 20 ml of 5.25 percent NaOCl for 15 minutes. It was delivered with the mas-



Fig. 1 MTA Angelus (Angelus, Londrina, Brazil) available in resealable vials. (Photos/Provided by Gary Glassman, DDS, FRCD(C))



Fig. 2 Radiograph of a necrotic lower left second premolar with large periradicular radiolucency with an incompletely formed root, both longitudinally and laterally.

tion.¹⁵ The development of normal, sterile granulation tissue within the root canal is thought to aid in revascularization and stimulation of cementoblasts or the undifferentiated mesenchymal cells at the periapex, leading to the deposition of a calcific material at the apex as well as on the lateral dentinal walls.¹²

tured during function, resulting in a microexposure and necrosis of the pulp. The tooth had an open apex associated with a large radiolucency (Fig. 2).

Periodontal probings were within normal limits for all teeth in the lower left region. Diagnostic testing was negative to cold and electric pulp

ter delivery tip and the macro canulae of the EndoVac apical negative pressure delivery system (Axis/SybronEndo, Coppel, Texas) (Fig. 3).

The canal was dried with paper points, and a mixture of ciprofloxacin, metronidazole and minocycline paste as described by Hoshino et al.¹⁷ was prepared into a creamy




Fig. 3 EndoVac apical negative pressure delivery system (Axis/ SybronEndo, Coppel, Texas)



Fig. 4 After the triple antibiotic paste was inserted into the canal, a temporary restoration was placed



Fig. 5 Blood clot was induced and MTA Angelus (Angelus, Londrina, Brazil) was placed over top, and then the tooth was restored with bonded composite



Fig. 6 Three-month recall reveals excellent longitudinal apical and lateral dentin development



Fig. 7 One-year recall radiograph reveals that definitive endodontics had been completed by the patient's new dentist

consistency and spun down the canal with a lentulo spiral instrument to a depth of 8 mm into the canal. The access cavity was closed with a sterile cotton pellet placed in the chamber and blue Cosmecore (Cosmedent, Chicago) (Fig. 4).

The patient returned three weeks later and was asymptomatic. The access was opened and the canal again flushed with 20 ml of 5.25 percent NaOCl for 15 minutes. It was delivered in the same manner as in the first visit with the master delivery tip and the macro canulae of the EndoVac apical negative pressure delivery system.

The canal appeared clean and dry, with no signs of inflammatory exudate. A #50 K-file was introduced into the canal until vital tissue was felt at a depth of 10 mm into the canal space. It was used to irritate the tissue gently to create some bleeding into the canal. The bleeding was stopped at a level of 5 mm below the level of the CEJ and left for 30 minutes, so that the blood would clot at that level.

After 30 minutes, the presence of the blood clot to approximately 5 mm apical of the CEJ was confirmed. White mineral trioxide aggregate, MTA Angelus was carefully placed over the blood clot and allowed to set for 20 minutes. After confirmation was achieved of its set, a bonded composite was placed and the patient was scheduled for follow-up in three months. Unfortunately, the MTA was placed further apically then would have been preferred (Fig. 5).

At the three-month follow-up appointment, the patient was totally asymptomatic, and the radiograph showed complete resolution of the radiolucency,

with closure of the apex and thickening of the dentinal walls. Pulp testing was inconclusive (Fig. 6).

At the one-year follow-up appointment, the radiograph revealed that treatment had been performed on this tooth by another dentist, different from her original dentist who made the initial referral. The new dentist, not familiar with revascularization treatment performed, had entered the root canal space, cleaned it out and obturated it with gutta-percha and sealer. Fortunately, the treatment was successful (Fig. 7).

Conclusion

The future of endodontics is bright as we continue to develop new techniques and technologies that will allow us to perform treatment painlessly and predictably and continue to satisfy one of the main objectives in dentistry — being to retain the natural dentition wherever possible and wherever practical.

References

1. Kakehashi S, Stanley HR, Fitzgerald RJ. The effects of surgical exposures of dental pulps in germ-free and conventional laboratory rats. *Oral Surg Oral Med Oral Pathol* 1965; 20: 340-349.
2. Moller AJR, Fabricius L, Dahlen G, Ohman A, Heyden G. Influence of periapical tissues of indigenous oral bacterial and necrotic pulp tissue in monkeys. *Scand J Dent Res* 1981; 89: 475-484.
3. Torabinejad M, Pitt Ford TR. Root end filling materials: a review. *Endod Dent Traumatol* 1996; 12:161-178.
4. Ribeiro DA. Do endodontic compounds induce genetic damage? A comprehensive review. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod*

2008; 105:251-256.

5. Enkel B, Dupas C, Armenogol V, et al. Bioactive materials in endodontics. *Expert Rev Med Devices* 2008; 5:475-494. that is hard tissue conductive (7).
6. Moretton TR, Brown CE Jr, Legan JJ, Kafrawy AH. Tissue reactions after subcutaneous and intraosseous implantation of mineral trioxide aggregate and ethoxybenzoic acid cement, hard tissue inductive, and biocompatible. *J Biomed Mater Res* 2000; 52:528-535.
7. Torabinejad M, Hong OU, Pitt Ford TR. Physical properties of a new root end filling material. *J Endodon* 1995; 21: 349-355.
8. Dentsply Tulsa Dental. Pro-Root™ MTA Root canal repair material; Material safety data sheet (MSDS).
9. Arnaldo Castellucci, MD, DDS. The Use of Mineral Trioxide Aggregate in Clinical and Surgical Endodontics. *Dentistry Today*, March 2003.

10. Duarte MA, Demarchi AC, Yamashita JC, Kuga MC, Fraga Sde C. pH and calcium ion release of 2 rootend filling materials. *Oral Surg Oral Med Oral Pathol Oral Radiol Endod.* 2003 Mar; 95(3):345-347.

Editorial note:

References are available from the author: DT

Contact Information

Gary Glassman, DDS, FRCD (C), graduated from the University of Toronto, Faculty of Dentistry in 1984; and graduated from the endodontology program at Temple University in 1987, where he received the Louis I. Grossman Study Club Award for academic and clinical proficiency in endodontics.

The author of numerous publications, Glassman lectures globally on endodontics, is on staff at the University of Toronto, Faculty of Dentistry, in the graduate department of endodontics, and is adjunct professor of dentistry and director of endodontic programming for the University of Technology, Jamaica.

He is a fellow of the Royal College of Dentists of Canada and the endodontic editor for *Oral Health dental journal*. He maintains a private practice, Endodontic Specialists, in Toronto, Ontario, Canada.

He can be reached through his website, www.rootcanals.ca.

mCME SELF INSTRUCTION PROGRAM

CAPP with Dental Tribune with its mCME- Self Instruction Program gives you the opportunity to have a quick and easy way to meet your continuing education needs. mCME offers you the flexibility to work at your own pace through the material from any location at any time. The content is international, drawn from the upper echelon of dental medicine, but also presents a regional outlook in terms of perspective and subject matter.

Membership:

Take membership for one year by subscription for the newspaper: 600 AED

Take article with one newspaper subscription: 100 AED per issue. After the payment, you will receive your membership number and will be able to start the program.

Completion of mCME

- mCME participants are required to read a continuing medical education (CME) article in each issue.
- Each article offers 2 CME Credit and followed by quiz questionnaire, which is available in <http://www.cappmea.com/mCME/questionnaires.html>.
- Each quiz has to be returned to events@cappmea.com or faxed to: +971436868883 in three months from the publication date
- A minimum passing score of 80% must be achieved in order to claim credit
- No more than two answered questions can be submitted in the same time
- Validity of the article – three months
- Validity of the subscription – one year
- Collection of Credit hours: you will receive the summary report with Certificate maximum one month after expire date of your membership. For single subscription Certificate and summary report will be send one month after the publication of the article.

The answers and critiques published herein have been checked carefully and represent authoritative opinions about the questions concerned.

Articles are available on www.cappmea.com after the publication. For more information please contact events@cappmea.com or +971 4 3616174

FOR INTERACTION WITH THE WRITERS FIND THE CONTACT DETAILS AT THE END OF EACH ARTICLE.

Earn 20 CME credit hours per year with annual membership in the mCME Program



Quick and easy CME points

- ✓ work at your own pace
- ✓ any location, any time

Accreditation
ADA CERP® Continuing Education Recognition Program
Continuing Education Recognition Program (CERP) is an ADA CERP Recognized Provider. ADA CERP is a service of the American Dental Association to assist dental professionals in identifying quality providers of continuing dental education. ADA CERP does not approve or endorse individual courses or instructors, nor does it imply endorsement of credit hours by boards of dentistry.

A new method for direct composite restoration of the posterior teeth

By Prof. Luca Giachetti, MD, DMD, MSc Department of Dentistry, University of Florence, Italy

Introduction

The evolution of composite materials and adhesive techniques has considerably changed the approach to restorations in posterior areas. The advantages of adhesive restorations are not only of an aesthetic nature, but, above all, relate to the possibilities of conserving a greater amount of healthy tissue and “reinforcing” the residual dental structure.

However, to exploit these advantages fully, we need rigorous clinical procedures which can limit what has always been the main flaw of composite materials: the polymerization shrinkage and the resulting stress which is responsible for most clinical failures.

Manufacturers have focused their efforts on producing materials which are ever easier to use and which, at the same time, are able to minimise their associated problems.

The recent introduction of the SonicFill™ System follows this direction. SonicFill combines the attributes of a low viscosity composite and a universal composite. By activating the composite with sonic energy, it is possible to fill the cavity and adapt the low viscosity material easily, and then compact and model it while the composite changes its consistency until it reaches a higher viscosity.

The manufacturer claims that it has the advantages of being: **Fast:** working time is reduced; it is possible to carry out single increments to an individual maximum thickness of 5 mm. **Reliable:** reduced shrinkage and good adaptability to the cavity walls due to the low initial viscosity.

Easy: it is possible to deliver the material using a small-diameter cannula and foot switch control. We present a clinical case below in which direct restorations have been produced with SonicFill on 3 elements of the 1° quadrant.

Clinical Case

Male patient, with an acceptable level of oral hygiene. In the maxillary right posterior quadrant, several deteriorated amalgam restorations are present with signs of marginal infiltration compatible with the age of the restorations, and signs of wear and tear in the zones of interocclusal contact. Tooth 1.5 has primary decay on the distal aspect of the tooth. The treatment plan was to replace the old amalgam restorations and to treat the primary caries with direct composites.



1. Initial case: 1.6 old amalgam with mesial-carries to be replaced, 1.5 primary distal decay, 1.4 old amalgam to be replaced



2. Isolation of the operative area with a rubberdam stabilised with a universal SoftClamp



3. Clinical situation after removal of the amalgam restorations. The contiguous elements are protected with metallic matrices before the marginal ridges are broken down



4. Access to the approximal carious lesions



5. Cavity cleaning, removal of demineralized tissue



6. Finishing of the margins with SonicSys inserts



7. Sectional metallic matrices contraposed on 1.6 and 1.5 stabilized with a wooden wedge, MetaFix All-in-One matrix stabilized with a wooden wedge on 1.4



8. Matrices in situ, the adaption at the level of the cervical floor can be seen



9. Application of the Self-Etch OptiBond XTR – Primer Adhesive System



10. Application of the Self-Etch OptiBond XTR – Bonding Adhesive System



11. Application of a thin layer of low viscosity Premise Flowable composite



12. Complete filling of the OM cavity of 1.6 with the SonicFill composite



13. Adaption of the material in the cavity with the CompoRoller oval tipped instrument and modelling with the point shaped tip




14. Application of the SonicFill composite in the occlusal distal cavity to the oblique ridge



15. Modelling the SonicFill composite with a Suter DD1-DD2 instrument



16. Restoration of the OD cavity of the 1.5 with SonicFill in a single application



17. Restoration of the OD cavity of the 1.4 with SonicFill in a single application





18. Characterization of the fissures with Kolor + Plus Brown



19. Restorations after finishing



20. Completed restorations after occlusion-check-up



21. Check-up at 1 week



15:39

15:42

15:47

15:52

15:55

Conclusion

The possibility of filling cavities to a depth of up to 5 mm with a single delivery effectively speeds up the work of performing composite restorations. The SonicFill composite presents good marginal adaption and is non-sticky. Once the sonic vibrations stop, it takes on an ideal consistency for modeling, and easily maintains the imposed shape. From an aesthetic point of view it is perhaps a little translucent to allow a greater depth of polymerization; however, it is possible to apply Kolor Plus®tints to make the restoration look natural. Ultimately, if the long-term controls show that the integrity of the margins is maintained, we will actually be able to confirm that we have accomplished a significant step towards simplifying direct restoration procedures with composite materials in posterior areas.

The products that appear in conjunction with this article are for illustrative or informational purposes only. Their inclusion does not denote endorsement by the author of this article. ^[1]

Contact Information



Prof. Luca Giachetti graduated in Medicine and Surgery in 1983 and specialized in Odontostomatology in 1986 at University of Florence Medical School. Msc in Dental Materials in 2009 at University of Siena Dental School. Chair of Dental Materials and Restorative Dentistry, University of Florence Dental School. Director of post graduate courses in aesthetics and adhesive dentistry, University of Florence Dental School. Dental Chief of Staff of Conservative Dentistry, Careggi Hospital-University, Italian NHS, Florence. He is member of the faculty in the International PhD program: "Biotechnology and Dental Biomaterials" at University of Siena Dental School. Associate Fellow in Education & Development, Warwick Dentistry, The University of Warwick, Coventry, UK. He has lectured at congresses and symposia and published on dental adhesives and composites in international dental journals. He runs a private practice in central Florence.

Just One Step for Bulk Filling A Giant Step for Dental Technology!

SonicFill is the only sonic-activated, single-step bulk filling system that enables you to go from placement to a polished restoration in less than 3 minutes. The change in viscosity allows perfect adaptation with optimal results and without loss of quality. Over 1,000,000 SonicFill tips have been sold and thousands of dentists are already using SonicFill every day in their practice. Kerr and KaVo innovating the way that dentist do posterior composite restorations!

Go to www.sonicfill.eu and see SonicFill in practice.
Or request a demonstration and join the revolution today!

SonicFill™



KerrHawe SA - Via Strecce 4 - CH-6934 Bioggio
Telephone +41 91 610 05 05 - www.kerrdental.eu

Your practice is our inspiration...

Kerr